

Dim. 125x39.5mm

Mr. INJECT™ DISPOSABLE SYRINGE
• PVC FREE • LATEX FREE • LUER MOUNT

5ml CE 2460

- Discard after single use - Sterilized by E.O.
- Do not use if previously opened or if pouch damaged.
- Do not use with Medicament that react with plastics.
- Contents Guaranteed Sterile, Non-Toxic, Non-Pyrogenic if unit pack is not opened or damaged.

PEEL OPEN

10°C 30°C

UDI

STERILE EO

Manufactured by:
BI-MED HEALTHCARE PRODUCTS PVT. LTD.
Unit-II, Plot No. 62-65, Sector -6A, SIDCUL, Haridwar 249 403, Uttarakhand (INDIA)
Mfg. Lic. No. MFG/MD/2020/000207
CE REP CMC Medical Devices & Drugs S.L. C/Horacio Lengua N18, CP29006, Málaga-Spain

REF : BM 0191
24GX1™
N. SIZE: (0.55x25mm)

LOT : 0522210
2022-07
2027-06

(01) 08906030481913

**100
Units**



**Mr. INJECT™
DISPOSABLE SYRINGE**

5ml

UDI

NEEDLE SIZE :



(01) 18906030481910(17) 000000(10) 000000(21) 000001

Expiry Date Batch Number Serial No. Box

LOT :



• NON-TOXIC • NON-PYROGENIC • STERILIZED BY EO • DISPOSABLE STERILE
• FOR SINGLE USE ONLY • STERILITY GUARANTEED TILL PACKAGE IS INTACT

Manufactured by:-

BIO-MED

HEALTHCARE PRODUCTS PVT. LTD.

Unit-II, Plot No. 62-65, Sector-6A, SIDCUL, Haridwar-249 403,

Uttarakhand (INDIA)

Mfg. Lic No. MFG/MD/2020/000207



UDI



STERILE EO



enquiry@biomedhealthcare.com

www.biomedhealthcare.com

ECIREP CMC Medical Devices & Drugs S.L. C/Horacio Lengo N18, CP29006, Málaga-Spain

DOCUMENT IDENTIFICATION NO. SY-02, ISSUE NO. 02, REV. NO. 02, DATE 23.07.2022

**1800
Units**




**Mr. INJECT™
DISPOSABLE SYRINGE**

5ml

UDI

NEEDLE SIZE :



(01) 58906030481918(17) 000000(10) 000000

Expiry Date Batch Number

LOT :



• NON-TOXIC • NON-PYROGENIC • STERILIZED BY EO • DISPOSABLE STERILE
• FOR SINGLE USE ONLY • STERILITY GUARANTEED TILL PACKAGE IS INTACT

Manufactured by:-

BIO-MED

HEALTHCARE PRODUCTS PVT. LTD.

Unit-II, Plot No. 62-65, Sector-6A, SIDCUL, Haridwar-249 403,

Uttarakhand (INDIA)

Mfg. Lic No. MFG/MD/2020/000207



UDI



STERILE EO



enquiry@biomedhealthcare.com

www.biomedhealthcare.com

ECIREP CMC Medical Devices & Drugs S.L. C/Horacio Lengo N18, CP29006, Málaga-Spain

DOCUMENT IDENTIFICATION NO. SY-03, ISSUE NO. 02, REV. NO. 02, DATE 23.07.2022

Duplex Box 205x135x190mm

 <p>Label sticker</p>	<p style="text-align: right;">100 Units</p> <p style="text-align: center;">DISPOSABLE SYRINGES</p> <p style="text-align: center;">5ml</p> <p style="text-align: center;">Safest Syringe</p> <ul style="list-style-type: none"> • PVC Free • Latex Free • Micro Fine Beveled Needle Tip • Smooth Surface Needle For Easy Gliding 	<p style="text-align: center;">DISPOSABLE SYRINGES</p> <ul style="list-style-type: none"> • STORE IN CLEAN, DRY AND INSECT FREE PLACE. • HANDLE WITH CARE. • BEFORE USE ROTATE THE NEEDLE CLOCKWISE TO ENSURE PROPER FITMENT ON THE BARREL NOZZLE. • PULL THE NEEDLE CAP STRAIGHT OFF TO AVOID DAMAGING NEEDLE POINT. • CHECK THE SYRINGE FOR VISUAL MANUFACTURING DEFECTS/CRACK BEFORE USING, TO AVOID LOSS OF MEDICATION. • AVOID EXPOSURE TO EXCESSIVE HEAT, ACIDIC/CORROSIVE FUMES. • PLEASE DISCARD IN SHARPS CONTAINER. <p>WARNING : • DO NOT USE IF PREVIOUSLY OPENED OR IF POUCH DAMAGED.</p> <p>STERILE EO</p> <p>STERILE • NON-TOXIC • NON-PYROGENIC • LUER MOUNT</p>	<p style="text-align: right;">100 Units</p> <p style="text-align: center;">DISPOSABLE SYRINGES</p> <p style="text-align: center;">5ml</p> <p style="text-align: center;">Safest Syringe</p> <ul style="list-style-type: none"> • PVC Free • Latex Free • Micro Fine Beveled Needle Tip • Smooth Surface Needle For Easy Gliding <p style="text-align: right; font-size: small;"> • THE PRODUCT SHOULD BE USED ACCORDING TO THE INSTRUCTIONS FOR USE. • FINISH BEFORE USE THIS PRODUCT / PACKING IS NOT THERMOSEAL OR DAMAGED. • FOR SINGLE USE ONLY. • STERILIZED BY EO. • HANDLE WITH CARE. • DISCARD AFTER USE. • DO NOT RECYCLE. • ENSURE PROPER DISPOSAL OF PRODUCT AND PACKING AFTER USE. • STORE IN COOL AND DRY PLACE. PROTECT FROM EXCESSIVE HEAT OR DIRECT SUNLIGHT. </p>
	<p style="text-align: center;">STERILE EO</p>		

<p style="text-align: center;">5ml</p> <p style="text-align: center;">DISPOSABLE SYRINGES</p> <p style="text-align: center;">100 Units</p>	<p style="text-align: center;">PRODUCTS AVAILABLE</p> <ul style="list-style-type: none"> • I.V. Cannula • I.V. Set • B.T. Set • 3-Way Stop Cock • 3-Way Stop Cock with Extension Tube • Extension Tube • Disposable Needle • Disposable Syringe • M.V. Set • Flow Regulator • Insulin Syringe
---	--

Carton Box_635x420x400mm


- FOR SINGLE USE ONLY.
- DISCARD AFTER USE. • DO NOT RECYCLE.
- STERILIZED BY EO.
- ENSURE PROPER DISPOSAL OF PRODUCT AND PACKING AFTER USE.
- THE PRODUCT SHOULD BE USED ACCORDING TO THE "INSTRUCTION FOR USE".
- ENSURE BEFORE USE THAT PRODUCT / PACKING IS NOT TAMPERED OR DAMAGED.

**1800
Units**

5ml

DISPOSABLE SYRINGES

**1800
Units**



5ml



Label sticker

serial no. of the Shipper/Carton

DISPOSABLE SYRINGES



• HANDLE WITH CARE • STORE IN A COOL AND DRY PLACE,
PROTECT FROM EXCESSIVE HEAT, DIRECT SUNLIGHT, RAIN OR EXCESSIVE HUMIDITY.

CARTON SIZE : 635x420x400 MM

- FOR SINGLE USE ONLY.
- DISCARD AFTER USE. • DO NOT RECYCLE.
- STERILIZED BY EO.
- ENSURE PROPER DISPOSAL OF PRODUCT AND PACKING AFTER USE.
- THE PRODUCT SHOULD BE USED ACCORDING TO THE "INSTRUCTION FOR USE".
- ENSURE BEFORE USE THAT PRODUCT / PACKING IS NOT TAMPERED OR DAMAGED.

**1800
Units**

5ml

DISPOSABLE SYRINGES

**1800
Units**

5ml

DISPOSABLE SYRINGES



• HANDLE WITH CARE • STORE IN A COOL AND DRY PLACE,
PROTECT FROM EXCESSIVE HEAT, DIRECT SUNLIGHT, RAIN OR EXCESSIVE HUMIDITY.

CARTON SIZE : 635x420x400 MM